

Basal Digital Thermometer



KD-233
(Flexible Tip)

KD-235

KD-269

Specifications

Model	KD-233	KD-235	KD-269
Measuring Range	89.60°F~109.39°F		
Resolution	0.01°C (0.01°F)		
Accuracy	±0.20°C (±0.40°F) as Temperature < 34.00°C (93.20°F) ±0.10°C (±0.20°F) as Temperature between 34.00°C ≤ T ≤ 42.00°C (93.20°F ≤ T ≤ 107.60°F) ±0.20°C (±0.40°F) as Temperature > 42.00°C (107.60°F)		
Display	Liquid Crystal Display 4 1/2 digits		
Operation Condition	Temperature: 10 to 40°C (50 to 104°F) Relative Humidity: 15 to 95 %		
Storage Condition	Temperature: -20 to 50°C (-4 to 120°F) Relative Humidity: 15 to 95 %		
Battery	Micro Alkaline 192 or LR 41 (1.55V)		
Power Consumption	0.15 milliwatt per measurement		
Low Battery Indication	Replace the battery if the low battery indicator appears		
Auto Off	Automatically after 10 minutes of non-use		
Size	126 x 19.9 x 9.3 mm	126.5 x 19.9 x 9.3 mm	133.8 x 17.3 x 9.4 mm
Weight	11g w/ battery	11.6g w/ battery	10.6g w/ battery
Standard	Complies with EN 12470-3 and ASTM E1112-00		

Specifications are subject to change for improvement without notice.

Accessories:

- Storage Case
- Probe Cover (optional)

K-jump Health Co., Ltd.

No. 56, Wu Kung 5th Rd., New Taipei Industrial Park,
New Taipei City 24890, Taiwan

Tel: 886-2-22991378 Fax: 886-2-22991386

E-mail: service@kjump.com.tw

<http://www.kjump.com.tw>



Features:



Regular



Waterproof



Memory



Beeper

Optional Functions:



Fever Alarm



C/F Switchable

