

K-jump®

Kinder Digital Thermometer (Rigid)



KD-121

Specifications

Model	KD-121
Measuring Range	32.0°C ~ 43.9°C(90.0°F~109.9°F)
Resolution	0.1°C (0.1°F)
Accuracy	±0.2°C(±0.4°F) as Temperature < 34.0°C(93.2°F) ±0.1°C(±0.2°F) as Temperature between 34.0°C ≤ T ≤ 42°C (93.2°F ≤ T ≤ 107.6°F) ±0.2°C(±0.4°F) as Temperature > 42.0°C(107.6°F)
Display	Liquid Crystal Display 3 1/2 digits
Operation Condition	Temperature: 10 to 40°C (50 to 104°F) Relative humidity: 15 to 95 %
Storage Condition	Temperature: -20 to 50°C (-4 to 120°F) Relative humidity : 15 to 95%
Battery	Micro Alkaline 192 or LR 41 (1.55V)
Power Consumption	0.15 milliwatt per measurement
Low Battery Indication	Replace the battery if the low battery indicator appears
Auto Off	Automatically after 9 minutes of non use
Size	154 x 21 x 12 mm
Weight	18 g w/ battery
Standard	Complies with EN 12470-3 and ASTM E1112-00

Specifications are subject to change for improvement without notice.

K-jump Health Co., Ltd.

No. 56, Wu Kung 5th Rd., New Taipei Industrial Park,
New Taipei City 24890, Taiwan

Tel: 886-2-22991378 Fax:886-2-22991386

E-mail: service@kjump.com.tw

<http://www.kjump.com.tw>



CE 0197

Features:



10" Fast Read



Beeper



Memory



Waterproof

Optional Functions:



Fever Alarm



C/F Switchable



HealthCare at Home
We Care for Your Health