



# Fast Read Digital Thermometer



## Features:

- 10" Fast Read
- Waterproof
- Beeper
- Memory
- Tip waterproof

## Optional Functions:

- Fever Alarm
- C/F Switchable

## Specifications

Model	KD-151	KD-152
Measuring Range	32.0°C ~ 43.9°C(90.0°F~109.9°F)	
Resolution	0.1°C (0.1°F)	
Accuracy	±0.2°C (±0.4°F) as Temperature < 34.0°C(93.2°F) ±0.1°C (±0.2°F) as Temperature between 34.0°C ≤ T ≤ 42°C (93.2°F ≤ T ≤ 107.6°F) ±0.2°C (±0.4°F) as Temperature > 42.0°C(107.6°F)	
Display	Liquid Crystal Display 3 1/2 digits	
Operation Condition	Temperature: 10 to 40°C ( 50 to 104°F) Relative humidity: 15 to 95 %	
Storage Condition	Temperature: -20 to 50°C(-4 to 120°F) Relative humidity : 15 to 95%	
Battery	Micro Alkaline 192 or LR 41 (1.55V)	
Power Consumption	0.15 milliwatt per measurement	
Low Battery indication	Replace the battery if the low battery indication appears	
Auto Off	Automatically after 9 minutes of non-use	
Size	140.5 x 17.8 x 10.2 mm	140 x 18 x 9 mm
Weight	12g	13g
Standard	Complies with EN 12470-3 and ASTM E1112-00	

Specifications are subject to change for improvement without notice.

## Accessories:

- Storage Case
- Probe Cover (optional)


## K-jump Health Co., Ltd.

No. 56, Wu Kung 5th Rd., New Taipei Industrial Park,  
New Taipei City 24890, Taiwan

Tel: 886-2-22991378 Fax:886-2-22991386

E-mail: [service@kjump.com.tw](mailto:service@kjump.com.tw)

<http://www.kjump.com.tw>



HealthCare at Home  
We Care for Your Health