


Fast Read Digital Thermometer



KD-181

Features:



Fast Read



Waterproof



Beeper



Memory



Flexible Tip

Optional Functions:



Fever Alarm



C/F Switchable

Specifications

Model	KD-181
Measuring Range	32.0°C ~ 43.9°C(90.0°F~109.9°F)
Resolution	0.1°C (0.1°F)
Accuracy	±0.2°C (±0.4°F) as Temperature < 34.0°C(93.2°F) ±0.1°C (±0.2°F) as Temperature between 34.0°C ≤ T ≤ 42°C (93.2°F ≤ T ≤ 107.6°F) ±0.2°C (±0.4°F) as Temperature > 42.0°C(107.6°F)
Display	Liquid Crystal Display 3 1/2 digits
Operation Condition	Temperature: 10 to 40°C (50 to 104°F) Relative humidity: 15 to 95 %
Storage Condition	Temperature: -20 to 50°C (-4 to 120°F) Relative humidity : 15 to 95%
Battery	Micro Alkaline 192 or LR 41 (1.55V)
Power Consumption	0.15 milliwatt per measurement
Low Battery indication	Replace the battery if the low battery indication appears
Auto Off	Automatically after 9 minutes of non-use
Size	129.5 x 21 x 12 mm
Weight	Approx. 12g
Standard	Complies with EN 12470-3 and ASTM E1112-00

Specifications are subject to change for improvement without notice.

Accessories:

- Storage Case
- Probe Cover (optional)

K-jump Health Co., Ltd.

No. 56, Wu Kung 5th Rd., New Taipei Industrial Park,
New Taipei City 24890, Taiwan

Tel: 886-2-22991378 Fax: 886-2-22991386

E-mail: service@kjump.com.tw

<http://www.kjump.com.tw>



CE 0197

