

K-jump®

Digital Pacifier Thermometer



KD-122



KD-122-1

Features:



Memory



Beeper



Waterproof

Optional Functions:



Fever Alarm



C/F Switchable

Specifications

| Model | KD-122 | KD-122-1 |
|------------------------|--|--------------------|
| Measuring Range | 32.0°C ~ 43.9°C(90.0°F~109.9°F) | |
| Resolution | 0.1°C (0.1°F) | |
| Accuracy | ±0.2°C (±0.4°F) as Temperature < 34.0°C(93.2°F) ±0.1°C (±0.2°F) as Temperature between 34.0°C ≤ T ≤ 42°C (93.2°F ≤ T ≤ 107.6°F) ±0.2°C (±0.4°F) as Temperature > 42.0°C(107.6°F) | |
| Display | Liquid Crystal Display 3 1/2 digits | |
| Operation Condition | Temperature: 10 to 40°C (50 to 104°F) Relative Humidity: 15 to 95 % | |
| Storage Condition | Temperature: -20 to 50°C (-4 to 120°F) Relative Humidity: 15 to 95 % | |
| Battery | Micro Alkaline 192 or LR 41 (1.55V) | |
| Power Consumption | 0.15 milliwatt per measurement | |
| Low Battery Indication | Replace the battery if the low battery indicator appears | |
| Auto Off | Automatically after 9 minutes of non-use | |
| Size | 58 x 38.4 x 46.6 mm | 53x 37.5 x 39.5 mm |
| Weight | Approx. 12g w/ battery | |
| Standard | Complies with EN 12470-3 and ASTM E1112-00 | |

Specifications are subject to change for improvement without notice.

Accessories:

- Storage Case

K-jump Health Co., Ltd.

No. 56, Wu Kung 5th Rd., New Taipei Industrial Park,
New Taipei City 24890, Taiwan


Tel: 886-2-22991378 Fax:886-2-22991386

E-mail: service@kjump.com.tw

<http://www.kjump.com.tw>



CE 0197



HealthCare at Home
We Care for Your Health