



Digital Pacifier Thermometer



KD-126

Features:



Beeper



Memory



Fever Alarm



Waterproof

Specifications

Model	KD-126
Measuring Range	32.0°C ~ 43.9°C(90.0°F~109.9°F)
Resolution	0.1°C (0.1°F)
Accuracy	±0.1°C for measuring range 35.5°C~42.0°C ±0.2°C for other measuring range
Display	Liquid Crystal Display 3 1/2 digits
Operation Condition	Temperature: 10 to 40°C (50 to 104°F) Relative Humidity: 15 to 95 %
Storage Condition	Temperature: -20 to 50°C (-4 to 120°F) Relative Humidity: 15 to 95 %
Battery	LR41 1.5V
Power Consumption	0.15 milliwatt per measurement
Low Battery indication	Replace the battery if the low battery indicator appears
Auto Off	Automatically after 9 minutes of non-use
Size	38.5 x 58.0 x 23.5 mm
Weight	Approx. 13g w/ battery
Standard	Complies with ASTM E1112-00

Specifications are subject to change for improvement without notice.

K-jump Health Co., Ltd.

No. 56, Wu Kung 5th Rd., New Taipei Industrial Park,
New Taipei City 24890, Taiwan


Tel: 886-2-22991378 Fax:886-2-22991386

E-mail: service@kjump.com.tw

<http://www.kjump.com.tw>



CE 0197



HealthCare at Home
We Care for Your Health