

Predictive Digital Thermometer



KD-195

- 4-second reading
- Jumbo LCD

Features:

-  Quick Read
-  Memory
-  Beeper
-  Fever Alarm
-  Three Mode Selection
-  Flexible Tip
-  Waterproof

Optional Functions:

-  C/F Switchable

Specifications

Model	KD-195
Measuring Range	30.5°C ~ 44.0°C(86.9°F~111.2°F)
Resolution	0.1°C (0.1°F)
Accuracy	±0.1°C (0.2°F) :34.0~42.0°C(93.2~107.6°F) ±0.2°C (0.4°F) :other range
Display	Liquid Crystal Display 3 1/2 digits
Operation Condition	Temperature: 10~40°C (50~104°F) Relative Humidity: 15~95 % non condensing
Storage Condition	Temperature: -20~60°C (-4~140°F) Relative Humidity: 15~95 % non condensing
Battery	CR 1632 (3.0V) x1
Power Consumption	6.0 milliwatt per measurement
Low Battery indication	Replace the battery if the low battery indicator appears
Auto Off	Automatically after 30-second of non-use
Size	134.8 mm x 23.8 mm x 12.1 mm
Weight	Approx. 16.2 grams including battery
Standard	Complies with EN 12470-3 and ASTM E 1112-00

Specifications are subject to change for improvement without notice.

Accessories:

- Storage Case
- Probe Cover (optional)

K-jump Health Co., Ltd.

No. 56, Wu Kung 5th Rd., New Taipei Industrial Park,
New Taipei City 24890, Taiwan

Tel: 886-2-22991378 Fax:886-2-22991386

E-mail: service@kjump.com.tw

<http://www.kjump.com.tw>



CE 0197

