

Flexible Tip Digital Thermometer



KD-133

Features:



Regular



Flexible Tip



Memory



Waterproof



Beeper

Optional Functions:



Fever Alarm



C/F Switchable

Specifications

Model	KD-133
Measuring Range	32.0°C ~ 43.9°C (90.0°F ~ 109.9°F)
Resolution	0.1°C (0.1°F)
Accuracy	±0.1°C (0.2°F): 34.0 ~ 42.0°C (93.2 ~ 107.6°F) ±0.2°C (0.4°F): other range
Display	Liquid Crystal Display 3 1/2 digits
Operation Condition	Temperature: 10 ~ 40°C (50 ~ 104°F) Relative Humidity: 15 ~ 95% non condensing
Storage Condition	Temperature: -20 ~ 60°C (-4 ~ 140°F) Relative Humidity: 15 ~ 95% non condensing
Battery	Micro Alkaline 192 or LR 41 (1.55V)
Power Consumption	Around 0.15 milliwatt per measurement
Low Battery Indication	Replace the battery if the low battery indicator appears
Auto Off	Automatically after 9 minutes of non-use
Size	127 x 20 x 10 mm
Weight	11 g w/ battery
Standard	Complies with ASTM-E1112, EN12470-3, EN60601-1

Specifications are subject to change for improvement without notice.

Accessories:

- Storage Case
- Probe Cover (optional)

K-jump Health Co., Ltd.

No. 56, Wu Kung 5th Rd., New Taipei Industrial Park,
New Taipei City 24890, Taiwan

Tel: 886-2-22991378 Fax: 886-2-22991386

E-mail: service@kjump.com.tw

<http://www.kjump.com.tw>



CE 0197

