

K-jump®

Digital Thermometer



KD-165

Features:



Optional Functions:



Specifications

Model	KD-165
Measuring Range	32.0°C ~ 43.9°C(90.0°F~109.9°F)
Resolution	0.1°C (0.1°F)
Accuracy	±0.1°C for measuring range 35.5°C~42.0°C ±0.2°C for other measuring range
Display	Liquid Crystal Display 3 1/2 digits
Operation Condition	Temperature: +10°C ~ +40°C (50°F ~ 104°F) Relative Humidity: 15% ~ 95%
Storage Condition	Temperature: -20 to 50°C (-4 to 120°F) Relative Humidity: 15 to 95 %
Battery	Micro Alkaline 192 or LR 41 (1.55V)
Power Consumption	0.15 milliwatt per measurement
Low Battery indication	Replace the battery if the low battery indicator appears
Auto Off	Automatically after 9 minutes of non-use
Size	146x21.8x11.2 mm
Weight	13.4g
Standard	Complies with EN-12470-3 and ASTM E1112-00

Specifications are subject to change for improvement without notice.

Accessories:

- Storage Case
- Probe Cover (optional)

K-jump Health Co., Ltd.

No. 56, Wu Kung 5th Rd., New Taipei Industrial Park,
New Taipei City 24890, Taiwan


Tel: 886-2-22991378 Fax:886-2-22991386

E-mail: service@kjump.com.tw

<http://www.kjump.com.tw>



CE 0197



HealthCare at Home
We Care for Your Health