

K-jump®

Compressor Nebulizer



KN-9330

Specifications

Power Supply	100~240V, 50~60Hz; 15Vdc, 1.2A
Power Consumption	Approx. 12W
Sound Level	Below 70dB
Particle Range(MMAD)	3.0~4.0 μ m
Ampoule Capacity	7ml
Flow Rate	0.3ml/min (NaCl 0.9%)
Intermittent Use	20 min ON/40 min OFF
Operation Conditions	+10C ~ +40C (50F~104 F); RH 10~95%
Storage Conditions	-25C ~ +60C (-13F ~ 140 F); RH 10~95%
Dimensions	Approx. 87 mm x 48 mm x 128.5 mm Approx. 87 mm x 48 mm x 186.5 mm(w/batt.)
Weight	Approx. 368g (w/o batt.) Approx. 610g (w/batt.)

Specifications are subject to change for improvement without notice.

Accessories:

- Pediatric Mask x 1
- Adult Mask x 1
- Mouthpiece x 1
- Carrying Bag x 1
- Air tubing x 1
- Connector x 1
- Mask Adapter x 1
- Inlet Filter x 5
- Outlet Filter x 5
- Rubber Cushion x 1
- Drug Ampoule Assembly x 1
- Rechargeable Battery Pack(optional) x 1
- Switching Adapter x 1
- Car Adapter (optional) x 1

K-jump Health Co., Ltd.

No.56 Wu Kung 5th Rd., WuKu Industrial Park,
Taipei Hsien, Taiwan 248
Tel: 886-2-22991378 Fax:886-2-22991386
E-mail: service@kjump.com.tw
<http://www.kjump.com.tw>



CE 0197



Features:

- Light weight and portable
- MMAD within respirable range
- Low battery indicator
- Safe & convenient



HealthCare at Home
We Care for Your Health