

# K-jump®

## Semiautomatic Upper Arm Blood Pressure Monitor



KP-6822



KP-6824



KP-6826



KP-6721

### Specifications

Model	KP-6721	KP-6822	KP-6824	KP-6826
Measuring Range	20 ~ 300 mmHg 40 ~ 200 pulse/minute			
Resolution	1 mmHg			
Accuracy	±3 mmHg or 2% (pressure) ±5% of readings (pulse rate)			
Measuring Method	Oscillometric method			
Operation Key	Two Keys: Power/Start, Memory			
Display	Liquid Crystal Display			
Inflation	Inflation bulb			
Exhaust	Manual release valve			
Memory	7-set, 42-set or 85-set			
Time	Displaying date, time and year			
Operation Condition	Temperature: 10°C ~ 40°C (50°F ~ 104°F) Relative Humidity: less than 85%			
Storage Condition	Temperature: -10°C ~ +50°C (+14°F ~ +122°F) Relative Humidity: less than 85%			
Power Source	6Vdc; "AA" type battery x 4			
Battery Life	Around 2000 measurements(alkaline batteries)			
Low Battery Indication	Replace the batteries if the low battery indicator appears			
Auto Off	Automatically after 150 seconds without operation			
Size	156 x 167 x 54 mm	140 x 80 x 47 mm	140 x 80 x 47 mm	140 x 80 x 47 mm
Weight	520g w/ batt.	244g w/ batt.	244g w/ batt.	244g w/ batt.
Standard	Complies with EN 1060-1, -3 and ANSI/AAMI SP10:2002			

Specifications are subject to change for improvement without notice.

### K-jump Health Co., Ltd.

No.56 Wu Kung 5<sup>th</sup> Rd., WuKu Industrial Park,  
Taipei Hsien, Taiwan 248

Tel: 886-2-22991378 Fax: 886-2-22991386

E-mail: [service@kjump.com.tw](mailto:service@kjump.com.tw)

<http://www.kjump.com.tw>



CE 0197

HealthCare at Home

We Care for Your Health

