

# K-jump®

## Automatic Upper Arm Blood Pressure Monitor



KP-6823

### Specifications

|                        |   |
|------------------------|---|
| Model                  | KP- 6823  |
| Measuring Range        | 20 ~ 300 mmHg<br>40 ~ 200 pulse/minute  |
| Resolution             | 1 mmHg  |
| Accuracy               | ±3 mmHg or 2% (pressure)<br>±5% of readings (pulse rate)                        |
| Measuring Method       | Oscillometric method  |
| Operation Key          | Two Keys: ON/OFF, Memory  |
| Display                | Liquid Crystal Display  |
| Inflation              | Micro rolling pump  |
| Exhaust                | Electrical Solenoid Valve   |
| Memory                 | 42-sets, or 85-sets   |
| Time                   | Displaying date, time and year  |
| Operation Condition    | Temperature: 10°C ~ 40°C (50°F ~ 104°F)<br>Relative Humidity: less than 85%     |
| Storage Condition      | Temperature: -10°C ~ +50°C (+14°F ~ +122°F)<br>Relative Humidity: less than 85% |
| Power Source           | 6Vdc; "AA" type battery x 4 or DC power input                                   |
| Battery Life           | Around 250 measurements (alkaline batteries)                                    |
| Low Battery Indication | Replace the batteries if the low battery indicator appears                      |
| Auto Off               | Automatically after 150 seconds without operation                               |
| Size                   | 140 x 80 x 46 mm  |
| Weight                 | Approx. 330g (including batteries)  |
| Standard               | Complies with EN 1060-1, -3 and ANSI/AAMI SP10:2002                             |

Specifications are subject to change for improvement without notice.

#### Accessories:

AC/DC adapter (optional)

### K-jump Health Co., Ltd.

No.56 Wu Kung 5<sup>th</sup> Rd., WuKu Industrial Park,  
Taipei Hsien, Taiwan 248



Tel: 886-2-22991378 Fax: 886-2-22991386

E-mail: [service@kjump.com.tw](mailto:service@kjump.com.tw)

<http://www.kjump.com.tw>



CE 0197



HealthCare at Home

We Care for Your Health