



# Wrist Watch Blood Pressure Monitor



KP- 6155



KP- 6165


## Specifications

Model	KP- 6155	KP- 6165
Measuring Range	20 ~ 300 mmHg 40 ~ 200 pulse/minute	
Resolution	1 mmHg	
Accuracy	±3 mmHg or 2% (pressure) ±5% of readings (pulse rate)	
Measuring Method	Oscillometric method	
Operation Key	Three Keys: Power/Start, Memory, Time Setting	
Display	Liquid Crystal Display	
Inflation	Micro rolling pump	
Exhaust	Electrical Solenoid Valve	
Memory	42-set, or 85-set	
Time	Displaying date, time and year	
Operation Condition	Temperature: 10°C ~ 40°C (50°F ~ 104°F) Relative Humidity: less than 85%	
Storage Condition	Temperature: -10°C ~ +50°C (+14°F ~ +122°F) Relative Humidity: less than 85%	
Power Source	3Vdc; "AAA" type battery x 2	
Battery Life	Around 250 measurements (alkaline batteries)	
Low Battery Indication	Replace the batteries if the low battery indicator appears	
Auto Off	Automatically after 150 seconds without operation	
Size	78 x 72 x 30 mm	78 x 72 x 30 mm
Weight	130g w/ batt.	
Standard	Complies with EN 1060-1, -3 and ANSI/AAMI SP10:2002	

Specifications are subject to change for improvement without notice.

## K-jump Health Co., Ltd.

No.56 Wu Kung 5<sup>th</sup> Rd., WuKu Industrial Park,  
 Taipei Hsien, Taiwan 248  
 Tel: 886-2-22991378 Fax:886-2-22991386  
 E-mail: service@kjump.com.tw  
 http://www.kjump.com.tw



HealthCare at Home  
 We Care for Your Health