

K-jump®

Talking Digital Thermometer



KV-171

Specifications

| | |
|------------------------|--|
| Model | KV-171 |
| Measuring Range | 32.0°C ~ 43.9°C(90.0°F~109.9°F) |
| Resolution | 0.1°C (0.1°F) |
| Accuracy | ±0.2°C(±0.4°F) as Temperature < 34.0°C(93.2°F) ±0.1°C(±0.2°F) as Temperature between 34.0°C ≤ T ≤ 42°C (93.2°F ≤ T ≤ 107.6°F) ±0.2°C(±0.4°F) as Temperature > 42.0°C(107.6°F) |
| Display | Liquid Crystal Display 3 1/2 digits |
| Operation Condition | Temperature: 10 to 40°C (50 to 104°F) Relative Humidity: 15 to 95 % |
| Storage Condition | Temperature: -20 to 50°C (-4 to 120°F) Relative Humidity: 15 to 95 % |
| Battery | Micro Alkaline 186 or LR 43 (1.55V) X 2 |
| Power Consumption | 0.4 milliwatt per measurement |
| Low Battery Indication | Replace the battery if the low battery indicator appears |
| Auto Off | Automatically after 10 minutes of non use |
| Size | 124 x 23.7 x 13 mm |
| Weight | Approx. 14g w/ battery |
| Standard | Complies with EN 12470-3 and ASTM E1112-00 |

Specifications are subject to change for improvement without notice.

Accessories:

- Storage Case
- Probe Cover (optional)

K-jump Health Co., Ltd.

No. 56, Wu Kung 5th Rd., New Taipei Industrial Park,
New Taipei City 24890, Taiwan

Tel: 886-2-22991378 Fax:886-2-22991386

E-mail: service@kjump.com.tw

<http://www.kjump.com.tw>




CE 0197

Features:



Optional Functions:



HealthCare at Home
We Care for Your Health